

Brizzy Patient Leaflet



This leaflet is intended to help the patient. For information regarding the limitations, contraindications, precautions, and warnings associated with Brizzy, please consult your doctor.

You can watch the video instructions here:



BRIZZY is a user-friendly portable sleep monitor. You can easily install the device yourself by following the instructions given here.


Installation of the Jawac sensor

Place each of the 2 probes on the velcro of an NS910 adhesive (JawaGrip) and fix the probes with the velcro:

1. Fasten the probe by closing the velcros.




2. Use one side of the velcro to secure the cable.




Both probes will be located on the same vertical axis and the cable must be placed on the same side. Make sure that the probes are secured to the face and cannot move or revolve.

Make sure that the wires are positioned so that they pass behind the same ear and adjust them by tightening the tubing around the cables.



Installation of the Device



Slip the abdominal belt into the slots of the device. Install the device on the belly (over a T-shirt or a pyjama shirt) by fastening the belt comfortably.

Starting - Up

Connect the distance sensor to the black socket of the device, the one marked «JAWAC».

Make sure you align the notch in the sensor connector with the crest on the socket before pushing in.



The device switches on upon connection of the sensors to it. A green indicator light starts flashing while the recording is in progress.

In the morning, upon waking up

Disconnect the connector of the black sensor (Jawac) and release the probes:

- Carefully loosen the velcro to remove the probes without pulling on the cables.
- Remove the adhesive from the face.

Take the equipment back to your doctor.



Nomics
Rue des Chasseurs
Ardennais 4
4031 Angleur (Belgium)
support@nomiscare.com

